



Product Datasheet

TINP1 Antibody (orb75404)

Description	TINP1 Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O95478
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...



Immunohistochemical analysis of formalin...