



Product Datasheet

PIAS3 Antibody (orb75419)

Description	PIAS3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y6X2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of K562 cell lysat...