


Product Datasheet

SIRT7 Antibody (orb75433)

Description	SIRT7 Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NRC8
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...