

Product Datasheet

Human Sulfite oxidase protein [Out of stock] (orb754715)

Description	Recombinant human Sulfite oxidase protein
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2um filtered solution in PBS with 5% trehalose, pH7.4
Storage	-20°C for 12 months as lyophilized;2-8°C for 1 month under sterile conditions after reconstitution
Tag	His tag
Note	For research use only
Purity	Greater than 95% by SDS-PAGE gel analyses.
Source	E.Coli
MW	51.2 kDa
Uniprot ID	P51687
Expiration Date	6 months from date of receipt.