



Product Datasheet

VPS39 Antibody (orb75569)

Description	VPS39 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96JC1
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat liver tissu...



Immunohistochemical
analysis of formalin...