


Product Datasheet

Human HN1 protein (orb756334)

Description	Recombinant human HN1 protein
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Storage	-20°C for 12 months as lyophilized. 2-8°C for 1 week under sterile conditions after reconstitution. Avoid repeated freeze-thaw cycles.
Tag	N-terminal His Tag
Note	For research use only
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Formula	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4
Purity	Greater than 95% as determined by SDS-PAGE
Source	E.Coli
MW	21 kDa
Uniprot ID	Q9UK76
Expiration Date	6 months from date of receipt.



Recombinant Human HN1 protein.