

Product Datasheet

PE2R4 antibody (orb756489)

Description	Rabbit polyclonal antibody to PE2R4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	PGE2 receptor EP4 subtype
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa, 22 kDa
Uniprot ID	P35408
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.