

Product Datasheet

GRM2 antibody (orb756514)

Description	Rabbit polyclonal antibody to GRM2
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	glutamate receptor, metabotropic 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	96 kDa
Uniprot ID	Q14416
Dilution Range	IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.