

Product Datasheet

NFKB1 antibody (orb756519)

Description	Mouse monoclonal antibody to NFKB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	105 kDa/50 kDa
Uniprot ID	P19838
Dilution Range	WB: 1:500-1:3000, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.