

Product Datasheet

TNFAIP3 antibody (orb756532)

Description	Rabbit polyclonal antibody to TNFAIP3
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	tumor necrosis factor, alpha-induced protein 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	90 kDa, 40 kDa
Uniprot ID	P21580
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.