

Product Datasheet

MYPC3 antibody (orb756569)

Description	Rabbit polyclonal antibody to MYPC3
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Myosin-binding protein C, cardiac-type
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	141 kDa
Uniprot ID	Q14896
Dilution Range	IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.