

Product Datasheet

MYD88 antibody (orb756632)

Description	Mouse monoclonal antibody to MYD88
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	myeloid differentiation primary response gene
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	35 kDa
Uniprot ID	Q99836
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.