

Product Datasheet

NeuN antibody (orb756633)

Description	Mouse monoclonal antibody to NeuN
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	hexaribonucleotide binding protein 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	A6NFN3
Dilution Range	IHC-P: 1:2000-1:20000
Expiration Date	12 months from date of receipt.