

Product Datasheet

IL1 beta antibody (orb756652)

Description	Rabbit polyclonal antibody to IL1 beta
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	interleukin 1, beta
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q63264
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.