

Product Datasheet

2019 nCOV S1 protein antibody (orb756658)

Description	Human polyclonal antibody to 2019 nCOV S1 protein
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	S1-RBD
Form/Appearance	liquid
Storage	PBS pH 7.4, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P0DTC2
Dilution Range	ELISA: 0.1-1 μg/ml, WB: 1:1000-1:10000
Expiration Date	12 months from date of receipt.