


Product Datasheet

Canine LH ELISA Kit (orb756666)

Description	Canine LH (Luteinizing Hormone) ELISA Kit
Reactivity	Canine
Range	0.781-50ng/ml
Target	LH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P18842
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Canine
LH(Luteinizing
Hormone)
ELISA Kit