


## Product Datasheet

### Human PML ELISA Kit (orb756686)

|                        |                                   |
|------------------------|-----------------------------------|
| <b>Description</b>     | Human PML (Protein PML) ELISA Kit |
| <b>Reactivity</b>      | Human                             |
| <b>Range</b>           | 78.125-5000pg/ml                  |
| <b>Target</b>          | PML                               |
| <b>Storage</b>         | 2-8 °C for 6 months               |
| <b>Note</b>            | For research use only             |
| <b>Uniprot ID</b>      | <b>P29590</b>                     |
| <b>Sensitivity</b>     | 46.875pg/ml                       |
| <b>Expiration Date</b> | Please enquire.                   |



Human  
PML(Protein  
PML) ELISA  
Kit