


## Product Datasheet

### Human PTPRS ELISA Kit (orb756689)

|                        |                                                                      |
|------------------------|----------------------------------------------------------------------|
| <b>Description</b>     | Human PTPRS (Receptor-type tyrosine-protein phosphatase S) ELISA Kit |
| <b>Reactivity</b>      | Human                                                                |
| <b>Range</b>           | 0.156-10ng/ml                                                        |
| <b>Target</b>          | PTPRS                                                                |
| <b>Storage</b>         | 2-8 °C for 6 months                                                  |
| <b>Note</b>            | For research use only                                                |
| <b>Uniprot ID</b>      | <b>Q13332</b>                                                        |
| <b>Sensitivity</b>     | 0.094ng/ml                                                           |
| <b>Expiration Date</b> | Please enquire.                                                      |



Human  
PTPRS(Receptor-  
type tyrosine-  
prote...