


## Product Datasheet

### Human RSPO3 ELISA Kit (orb756695)

|                        |                                     |
|------------------------|-------------------------------------|
| <b>Description</b>     | Human RSPO3 (R-spondin-3) ELISA Kit |
| <b>Reactivity</b>      | Human                               |
| <b>Range</b>           | 15.625-1000pg/ml                    |
| <b>Target</b>          | RSPO3                               |
| <b>Storage</b>         | 2-8 °C for 6 months                 |
| <b>Note</b>            | For research use only               |
| <b>Uniprot ID</b>      | <b>Q9BXY4</b>                       |
| <b>Sensitivity</b>     | 9.375pg/ml                          |
| <b>Expiration Date</b> | Please enquire.                     |



Human  
RSPO3(R-  
spondin-3)  
ELISA Kit