


## Product Datasheet

### Human SARM1 ELISA Kit (orb756697)

|                        |                                                                          |
|------------------------|--------------------------------------------------------------------------|
| <b>Description</b>     | Human SARM1 (Sterile alpha and TIR motif-containing protein 1) ELISA Kit |
| <b>Reactivity</b>      | Human                                                                    |
| <b>Range</b>           | 0.156-10ng/ml                                                            |
| <b>Target</b>          | SARM1                                                                    |
| <b>Storage</b>         | 2-8 °C for 6 months                                                      |
| <b>Note</b>            | For research use only                                                    |
| <b>Uniprot ID</b>      | <b>Q6SZW1</b>                                                            |
| <b>Sensitivity</b>     | 0.094ng/ml                                                               |
| <b>Expiration Date</b> | Please enquire.                                                          |



Standard  
Curve of  
Human  
SARM1  
ELISA Kit.