


Product Datasheet

Human SPTBN1 ELISA Kit (orb756702)

Description	Human SPTBN1 (Spectrin beta chain, non-erythrocytic 1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SPTBN1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q01082
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
SPTBN1(Spectrin
beta chain, non-
er...