


Product Datasheet

Human SYNGAP1 ELISA Kit (orb756705)

Description	Human SYNGAP1 (Synaptic Ras GTPase-activating protein 1) ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	SYNGAP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96PV0
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
SYNGAP1(Synaptic
Ras GTPase-
activa...