


Product Datasheet

Human TNFAIP1 ELISA Kit (orb756709)

Description	Human TNFAIP1 (Tumor necrosis factor, alpha-induced protein 1, endothelial) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	TNFAIP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13829
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
TNFAIP1(Tumor
necrosis factor,
alp...