


Product Datasheet

Human ITIH3 ELISA Kit (orb756732)

Description	Human ITIH3 (Inter-alpha-trypsin inhibitor heavy chain H3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	ITIH3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q06033
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
ITIH3(Inter-
alpha-
trypsin
inhibito...