


Product Datasheet

Human PCBD2 ELISA Kit (orb756745)

Description	Human PCBD2 (Pterin-4-alpha-carbinolamine dehydratase 2) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	PCBD2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H0N5
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
PCBD2(Pterin-4-
alpha-
carbinolamine...