


Product Datasheet

Human Pdha1 ELISA Kit (orb756752)

Description	Human Pdha1 (Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial) ELISA K...
Reactivity	Human
Range	0.313-20ng/ml
Target	Pdha1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P08559
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
Pdha1(Pyruvate
dehydrogenase
E1 co...