


Product Datasheet

Human anti-SARS-CoV2 (S) IgG ELISA Kit (orb756760)

Description	Human anti-SARS-CoV2 (S) IgG ELISA Kit
Reactivity	Human
Range	3.906-250ng/ml
Target	anti-SARS-CoV2(S) IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	2.344ng/ml
Expiration Date	Please enquire.



Human
anti-
SARS-
CoV2(S)
IgG ELISA
Kit