


## Product Datasheet

### Human CDCP1 ELISA Kit (orb756767)

|                        |                                                            |
|------------------------|------------------------------------------------------------|
| <b>Description</b>     | Human CDCP1 (CUB domain-containing protein 1)<br>ELISA Kit |
| <b>Reactivity</b>      | Human                                                      |
| <b>Range</b>           | 0.156-10ng/ml                                              |
| <b>Target</b>          | CDCP1                                                      |
| <b>Storage</b>         | 2-8 °C for 6 months                                        |
| <b>Note</b>            | For research use only                                      |
| <b>Uniprot ID</b>      | <b>Q9H5V8</b>                                              |
| <b>Sensitivity</b>     | 0.094ng/ml                                                 |
| <b>Expiration Date</b> | Please enquire.                                            |



Human  
CDCP1(CUB  
domain-  
containing  
protei...