


## Product Datasheet

### Human CXorf61 ELISA Kit (orb756770)

|                        |                                                                |
|------------------------|----------------------------------------------------------------|
| <b>Description</b>     | Human CXorf61 (KITa-kyushu lung cancer antigen 1)<br>ELISA Kit |
| <b>Reactivity</b>      | Human                                                          |
| <b>Range</b>           | 78.125-5000pg/ml                                               |
| <b>Target</b>          | CXorf61                                                        |
| <b>Storage</b>         | 2-8 °C for 6 months                                            |
| <b>Note</b>            | For research use only                                          |
| <b>Uniprot ID</b>      | <b>Q5H943</b>                                                  |
| <b>Sensitivity</b>     | 46.875pg/ml                                                    |
| <b>Expiration Date</b> | Please enquire.                                                |



Human  
CXorf61(KITa-  
kyushu lung  
cancer an...