


## Product Datasheet

### Human ELAVL2 ELISA Kit (orb756773)

|                        |                                              |
|------------------------|----------------------------------------------|
| <b>Description</b>     | Human ELAVL2 (ELAV-like protein 2) ELISA Kit |
| <b>Reactivity</b>      | Human                                        |
| <b>Range</b>           | 0.156-10ng/ml                                |
| <b>Target</b>          | ELAVL2                                       |
| <b>Storage</b>         | 2-8 °C for 6 months                          |
| <b>Note</b>            | For research use only                        |
| <b>Uniprot ID</b>      | <b>Q12926</b>                                |
| <b>Sensitivity</b>     | 0.094ng/ml                                   |
| <b>Expiration Date</b> | Please enquire.                              |



Human  
ELAVL2(ELAV-  
like protein 2)  
ELISA ...