


## Product Datasheet

### Human GPR43 ELISA Kit (orb756784)

|                        |                                                    |
|------------------------|----------------------------------------------------|
| <b>Description</b>     | Human GPR43 (Free fatty acid receptor 2) ELISA Kit |
| <b>Reactivity</b>      | Human                                              |
| <b>Range</b>           | 0.156-10ng/ml                                      |
| <b>Target</b>          | GPR43                                              |
| <b>Storage</b>         | 2-8 °C for 6 months                                |
| <b>Note</b>            | For research use only                              |
| <b>Uniprot ID</b>      | <b>O15552</b>                                      |
| <b>Sensitivity</b>     | 0.094ng/ml                                         |
| <b>Expiration Date</b> | Please enquire.                                    |



Human  
GPR43(Free  
fatty acid  
receptor 2)

...