


Product Datasheet

Human KIR3DL1 ELISA Kit (orb756790)

Description	Human KIR3DL1 (Killer cell immunoglobulin-like receptor 3DL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	KIR3DL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P43629
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
KIR3DL1(Killer
cell
immunoglobulin...