


Product Datasheet

Mouse Anti-2019 nCoV (S) IgG ELISA Kit (orb756799)

Description	Mouse Anti-2019 nCoV (S) IgG ELISA Kit
Reactivity	Mouse
Range	1.563-100ng/ml
Target	Anti-2019 nCoV(S)IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Standard
Curve of
Mouse
Anti-2019
nCoV (...)