


Product Datasheet

Human AZGP ELISA Kit (orb756812)

Description	Rat AZGP1 (Zinc-alpha-2-glycoprotein) ELISA Kit
Reactivity	Rat
Range	3.13-200ng/ml
Target	AZGP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q63678
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Rat
AZGP1(Zinc-
alpha-2-
glycoprotein)
ELI...