


Product Datasheet

Rat beta Arrestin-1 ELISA Kit (orb756813)

Description	Rat beta Arrestin-1 (Beta Arrestin- 1) ELISA kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	β Arrestin-1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P29066
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat β
Arrestin-1
(Beta
Arrestin-1)
ELIS...