


Product Datasheet

Human beta Arrestin-2 ELISA Kit (orb756814)

Description	Rat beta Arrestin-2 (Beta Arrestin-2) ELISA kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	β Arrestin-2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P29067
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat β
Arrestin-2
(Beta
Arrestin-2)
ELIS...