


Product Datasheet

12-HETE ELISA Kit (orb756821)

Description	12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
Reactivity	All
Range	0.313-20ng/ml
Target	12-HETE
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



12-HETE (12-Hydroxyeicosatetraenoic Acid...)