

Product Datasheet

CDKN2A Antibody (orb757204)

Description	Mouse monoclonal antibody to CDKN2A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Purified recombinant fragment of human CDKN2A (AA: 1-156) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	16.5kDa
Entrez	1029
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200-1/1000; RCM: 1/200-1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.