

Product Datasheet

MYCN Antibody (orb757227)

Description	Mouse monoclonal antibody to MYCN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Purified recombinant fragment of human MYCN (AA: 337-464) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	49.5kDa
Entrez	4613
Dilution Range	IHC: 1/200-1/1000; RCM: 1/200-1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.