

Product Datasheet

SARS-Cov2-NP3 (orb757280)

Description	Mouse monoclonal antibody to SARS-Cov2-NP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of SARS-CoV-2-N (AA: 240-419) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	23.7kDa
Entrez	43740575
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.