



Product Datasheet

FLI1 Antibody (orb757315)


Description	Mouse monoclonal antibody to FLI1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Purified recombinant fragment of human FLI1 (AA: 303-452) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	50.9kDa
Entrez	2313
Dilution Range	IHC: 1/200-1/1000; RCM: 1/200-1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Hela cells u...



Line graph
illustrates
about the Ag-
Ab r...



Western blot
analysis of
recombinant
pro...