

Product Datasheet

FEN antibody [Out of stock] (orb757551)

Description	Mouse monoclonal antibody to Fentanyl
Species/Host	Mouse
Conjugation	Unconjugated
Immunogen	Fentanyl-KLH
Preservatives	PBS, pH 7.4*=>7.2
Concentration	4.8
Storage	Storage before 8 weeks: 2°C-8°C Storage until Expiration: -20°C +/- 10°C
Note	For research use only
Purity	> 95%
Expiration Date	12 months from date of receipt.