


## Product Datasheet

### DMG-PEG 2000 (orb757605)

<b>Description</b>	DMG-PEG 2000 is used for the preparation of liposome for siRNA delivery with improved transfection e...
<b>Note</b>	For research use only
<b>Formula</b>	$(C_2H_4O)_n C_{32}H_{62}O_5$
<b>CAS Number</b>	[160743-62-4]
<b>SMILES</b>	<chem>CCCCCCCCCCCC(OC(COCCOC)COC(CCCCCCCCCCCC)=O)=O.[n]</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical  
Structure  
(160743-62-4)