



Product Datasheet

CLEC2 antibody (orb757706)

Description	Mouse monoclonal antibody to CLEC2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP
Immunogen	A recombinant extracellular domain of human CLEC2 (amino acids 68-229)
Target	CLEC2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9P126
Entrez	51266
Dilution Range	FC: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human CLEC2 positive throm...



Flow cytometry surface staining pattern ...