



Product Datasheet

TCR Vdelta2 antibody (orb757707)


Description	Mouse monoclonal antibody to TCR Vdelta2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	human recombinant TCR Vdelta2
Target	TCR Vdelta2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	FC: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation
of human
TCR Vdelta2
positive...



Flow
cytometry
surface
staining
pattern ...