

Product Datasheet

TrkB (Phospho-Tyr705) antibody (orb758123)

Description	Rabbit polyclonal antibody to TrkB (Phospho-Tyr705)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q16620
Expiration Date	12 months from date of receipt.