

Product Datasheet

Calnexin (Phospho-Ser583) antibody (orb758135)

Description	Rabbit polyclonal antibody to Calnexin (Phospho-Ser583)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27824
Expiration Date	12 months from date of receipt.