

Product Datasheet

n-NOS (Phospho-Ser852) antibody (orb758144)

Description	Rabbit polyclonal antibody to n-NOS (Phospho-Ser852)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P29475
Expiration Date	12 months from date of receipt.