

Product Datasheet

14-3-3 Theta/Tau (Phospho-Ser232) antibody (orb758153)

Description	Rabbit polyclonal antibody to 14-3-3 Theta/Tau (Phospho-Ser232)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27348
Expiration Date	12 months from date of receipt.