

Product Datasheet

CD18/ITGB2 (Phospho-Thr758) antibody (orb758169)

Description	Rabbit polyclonal antibody to CD18/ITGB2 (Phospho-Thr758)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05107
Expiration Date	12 months from date of receipt.